






Omnican® fine

PEN NEEDLES FOR INSULIN INJECTION



Omnican® fine		4 mm	6 mm	8 mm	10 mm	12 mm
Dimensions		0.25 mm x 4 mm 31G x 5/32"	0.25 mm x 6 mm 31G x 1/4"	0.25 mm x 8 mm 31G x 5/16"	0.30 mm x 10 mm 30G x 3/8"	0.33 mm x 12 mm 29G x 1/2"
						
REF	Needle shield – yes	9167811EX	9168956EX	9168978EX	–	9168990EX
	Needle shield – no	9167820WE	9167668WE	9167684WE	9167625WE	9167641WE
Material	Needle tube	Stainless steel, silicone oil coated				
	Needle hub	Copolyester				
	Needle container	Polyethylene				
	Needle shield	Styrene-butadiene copolymer				
Wall thickness cannula		6 mm - 12 mm: Thin-wall according to ISO 9626:2016				
		4 mm: Extra-thin-wall according to ISO 9626:2016				
Sales unit		100 Omnican® fine pen needles in a dispenser box				
Shelf life		5 years (60 months)				
Sterilization		Ethylene oxide (EO)				
Certificates		EN ISO 13485:2012				
Product standards		EN ISO 11608-2:2012				
Product classification		Medical Device according to MDD 93/42/EEC, Class IIa				

WE = Europe | EX = Asia Pacific


























All data refer to typical single values and may be subject to alterations.

Omnican® fine

COMPATIBILITY CHART



Omnican® fine pen needles fit to major brands of insulin injectors according to ISO 11608-2:2012

Pen	Type	Manufacturer	Product Name	Pen	Type	Manufacturer	Product Name
	R	Berlin-Chemie	BerliPen® areo 2		R	Novo Nordisk	NovoPen Echo®
	R	Berlin-Chemie	BerliPen® Junior		R	Novo Nordisk	NovoPen® Junior
	R	Berlin-Chemie	BerliPen® Precision		D	Novo Nordisk	Victoza®
	D	BMS* AstraZeneca	Byetta®		R	Owen Mumford	Autopen® 24
	R	Lilly	HumaPen® Luxura		R	Owen Mumford	Autopen® Classic 1
	R	Lilly	HumaPen® Luxura HD		R	Owen Mumford	Autopen® Classic 2
	R	Lilly	HumaPen® Memoir		R	Sanofi	ClikStar®
	R	Lilly	HumaPen® Savvio		D	Sanofi	Lyxumia®
	D	Lilly	KwikPen™		D	Sanofi	SoloStar®
	D	Novo Nordisk	FlexPen®		R	Sanofi	TactiPen®
	D	Novo Nordisk	FlexTouch®		R	Ypsomed	ServoPen®
	R	Novo Nordisk	NovoPen® 4		R	Ypsomed	YpsoPen®
	R	Novo Nordisk	NovoPen® 5				

* BMS = Bristol-Myers Squibb | D = Disposable | R = Reusable | The brands listed are the property of the named companies.

All data refer to typical single values and may be subject to alterations.